

RETURN FORM M3400

V20240321

ITEMS WITH MACHINE

	R1 SLEEVE (MACH3500055)	1	CARRY CASE (CASE3400001)	
	B1 BOB (MACH3500058)		USB TO SERIAL ADAPTER (USBC0000001)	
	B2 BOB (MACH3500061)		SERIAL/DATA CABLE (CABL3500001)	
	DUST SHIELD (MACH3500057)		POWER SUPPLY (POWR3400001)	
	SAMPLE CUP 180 mL (CUP0000002)		THERMOMETER (THEM3400001)	
	HEATER CUP (HEAT3410AS1) (HEAT3410AS2)			
DATE IN:		MA	AACHINE ID #:	
DATE IN.		IVIA	VIACITINE ID #.	
	CLIENT I	INFOR	MATION	
CONTACT NAME:		PH	PHONE #:	
COMPANY NAME:				
EMAIL ADDRESS:				
COMPANY ADDRESS:				
REASON FOR RETURN				

When returning machines for service, please include the following items:

- 1. The B1 bob, but not installed on the shaft (shipping damage may occur).
- 2. The Dust Shield.
- 3. The Sleeve/Rotor.
- 4. The Temperature Probe.
- 5. Any accessories you want to have inspected/service (heater cup, additional bobs, sample cup).

NOTE: Pack the M3600 hardware carefully for shipment. Be sure that all included components are protected from damage. Grace Instrument is not responsible for damage to items during transit to our facility.